Quantitative determination of anti-streptolysin O (ASO) IVD

Store 2-8°C.

PRINCIPLE OF THE METHOD

The ASO-Turbilatex is a quantitative turbidimetric test for the measurement of ASO in human serum or plasma.

Latex particles coated with streptolysin O (SLO) are agglutinated when mixed with samples containing ASO. The agglutination causes an absorbance change, dependent upon the ASO contents of the patient sample that can be quantified by comparison from a calibrator of known ASO concentration.

CLINICAL SIGNIFICANCE

SLO is a toxic immunogenic exoenzyme produced by β -hemolytic Streptococci of groups A, C and G. Measuring the ASO antibodies are useful for the diagnostic of rheumatoid fever, acute glomerulonephritis and streptococcal infections. Rheumatic fever is an inflammatory disease affecting connective tissue from several parts of human body (skin, heart, joints, etc) and acute glomerulonephritis is a renal infection that affects mainly to renal glomerulus.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, pH 8,2. Preservative
Latex (R2)	Latex particles coated with streptolysin O, pH 10,0. Preservative.
ASO-CAL	Calibrator. Human serum. ASO concentration is stated on the vial label.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Use ASO Calibrator Reference 1107012.

The sensitivity of the assay and the target value of the calibrator have been standardized against the ASO International Calibrator (WHO).

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use.

ASO Calibrator: Reconstitute (→) with 1,0 mL of distilled water. Mix gently and incubate at room temperature for 10 minutes before use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations prevented during their use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles and turbidity.

ASO Calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Normal values up to 200 IU/mL (adults) and 100 IU/mL (children < 5 years old)⁶.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Controls ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240
Test	ASO	Vol. R2	60
Full Name	Anti-streptolysin O	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	UI/mL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Reagent Blank	20-21
Sec. Wave.		React. Time	44-60

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	ASO	No.	**
Full name	ASO	Print name	ASO
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	UI/mL	Decimal	0.1
Reagent Blank	20 -21	React. Time	54 - 70
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until 31 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 800 IU/mL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/3 in NaCl 9 g/L and retested again. The linearity limit depends on the sample-reagent ratio, as well the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 20 IU/mL give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 1000 IU/mL.
- Sensitivity:** Δ 0,73 mA. IU/mL.
- Precision:** The reagent has been tested for 20 days, using three different ASO concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 100 IU/mL	+/- 200 IU/mL	+/- 400 IU/mL
Total	6,4%	5,7%	5,1%
Within Run	2,4%	1,7%	1,4%
Between Run	3,6%	4,2%	4,9%
Between Day	4,7%	3,5%	0,7%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 60 samples of different concentrations of ASO were assayed. The correlation coefficient (r) was 0.99 and the regression equation $y = 0,915x - 4,844$. The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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APPLICATION AU SPIN640

TEST INFORMATION		REAGENT VOLUME	
N°	**	Vol. R1	240
Test	ASO	Vol. R2	60
Full Name	Anti-streptolysin O	Vol. R3	
Standard n°	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	UI/mL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Reagent Blank	20-21
Sec. Wave.		React. Time	44-60

APPLICATION AU SPIN640Plus

EDIT PARAMETERS			
Test	ASO	No.	**
Full name	ASO	Print name	ASO
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	UI/mL	Decimal	0.1
Reagent Blank	20-21	React. Time	54-70
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

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